

Murine Anti-bovine Factor X

Clone GMA-560

Bovine Factor X (Mr 55,000) is a vitamin K-dependent plasma protein zymogen that plays a central role as the substrate for both the intrinsic (factor VIIa, tissue factor) and extrinsic (factor IXa, factor VIIIa) pathways. In the presence of cofactor factor Va, phospholipid, and Ca^{2+} , activated factor X cleaves two peptide bonds in prothrombin to form thrombin. GMA-560 (also known as α -BFX-2b) binds bovine FX and factor Xa in solid-phase ELISA and Western blot. GMA-560 is partially calcium dependent, inhibits FX activity in human, bovine, porcine, rabbit and canine plasma, inhibits conversion of prothrombin to thrombin, and prevents inactivation of FXa by antithrombin III¹. It can be used in purification of factor X from plasma^{2,3} and blocks FXa-thrombomodulin interaction⁴.

Description

Antibody Source: mouse monoclonal, IgG₁

Antigen Species Bound: human, bovine, canine, rabbit, sheep, porcine¹

Specificity: factor X/Xa heavy chain

Immunogen: bovine factor X

Formulation and Storage

Purity: Purified by protein G affinity chromatography from serum-free cell culture supernatant.

Product Formulation: Lyophilized from a ≥ 1 mg/ml solution in 20 mM NaH_2PO_4 0.15 M NaCl, 1.0% (w/v) mannitol, pH 7.4. Concentration determined by absorbance measurement at 280 nm and using an extinction coefficient of 1.4 ($\epsilon_{0.1\%}$).

Reconstitution: Reconstitute with deionized water.

Storage: Store lyophilized or reconstituted and aliquoted material at -20°C for prolonged periods. Avoid freeze-thaw cycles. Alternatively, add 0.02% (w/v) sodium azide to reconstituted solution and store at 4°C .

Country of Origin: USA


Size Options: 0.1 mg or 0.5 mg

Applications

Working Concentration: Approximately 1-5 $\mu\text{g/ml}$. Researcher should titer antibody in specific assay.

ELISA: Binds bovine and human factor X and Xa in presence of calcium.

Immunoblotting: Binds bovine factor X and factor Xa heavy chain under reduced and non-reduced conditions in presence of calcium.



References

- [1] W.R. Church, T.L. Messier, M.M. Tucker, K. Mann. An inhibitory monoclonal antibody to factor X that blocks prothrombin activation but not prothrombinase enzyme assembly. (1988). *Blood*. 72:1911-1921.
- [2] W.R. Church, K. Mann. A simple purification of human factor X using a high affinity monoclonal antibody immunoabsorbant. (1985). *Thromb Res*. 38(4):417-424.
- [3] R. Jenny, W. Church, B. Odegaard, R. Litwiller, K. Mann. Purification of six human vitamin K-dependent proteins in a single chromatographic step using immunoaffinity columns. (1986). *Prep Biochem*. 16(3):227-245.
- [4] P.E. Haley, M.F. Doyle, K.G. Mann. The Activation of Bovine Protein C by Factor Xa.(1989). *J Biol Chem*. 264(27):16303-16310.
- [5] M. Wilkens, S. Krishnaswamy. The Contribution of Factor Xa to Exosite-dependent Substrate Recognition by Prothrombinase. (2002). *J Biol Chem*. 277(11):9366-9374.